

Am I Pregnant?™

QUESTIONS AND ANSWERS

1. **How does the test work?**

Am I Pregnant Pregnancy Test detects hCG (human Chorionoc Gonadotropin), a hormone in the urine once you become pregnant. A positive result indicates the hCG has been detected.

2. **What do I do if the test is positive?**

It means that the pregnancy hormone hCG has been found in your urine. You can assume you are pregnant and should see your doctor. It is recommended that you retest three days later to confirm.

3. **What do I do if the test is negative?**

It means that no pregnancy hormone has been detected and you are probably not pregnant. If your period does not start within a week, retest by using another Am I Pregnant Pregnancy Test. If the test is still negative, you should see your doctor.

4. **What do I do if the test is invalid?**

Improper operation or use of an expired test stick is the most frequent causes of an invalid result. In the event of an invalid result, it is recommended you repeat the test to by using a new test stick and fresh urine sample.

5. **Do I have to test in the morning?**

First morning urine usually contains a higher level of hCG. If you perform the test before expected period, it is recommended you use the first morning urine. Otherwise, you can test at any time of the day. Drinking excessive fluid before testing could dilute hCG in your urine and cause an inaccurate result.

6. **What if the test is allowed to run for more than 5 minutes?**

A positive result will not change except that the background may get darker. A negative result may show a faint band in T area. But, this is not a positive result. Do not interpret a result after 10 minutes.

7. **How accurate is the Am I Pregnant Pregnancy Test?**

In both laboratory and clinical studies with consumers, this test shown to be over 99% accurate* when used from the first day of your expected period.

8. **Can any medicines affect the result?**

Only medicines which contain hCG can give a false reading. HCG-containing drugs include: A.P.L., Pregnyl, Follutein and Profasi. Hormone therapies (such as Clomid, Serophene or Pergonal), painkillers, or the

contraceptive pill should not affect the test result. Drinking alcohol will not affect the test result.

9. **Can I get a false or misleading result?**

If you follow the procedure incorrectly, you may get false results. Certain rare medical conditions can also cause a false result. If you repeatedly get unexpected results, you should check with your doctor.

10. **How soon can I use Am I Pregnant Pregnancy Test?**

This test can detect hCG pregnancy hormone in your urine as early as 4 days before your expected period as well as any days thereafter. If you test before your expected period and get a negative, there is still a possibility you may be pregnant.

11. **I have an irregular period. How do I know I am testing at the right time?**

You should test based on your longest cycle in recent months. If you are not sure when your period is due we suggest testing at least 19 days from the day you think you became pregnant. If you still read negative you should retest in 3 days because you may have tested too early.

12. **What if I splash the result window with urine?**

If the window and the reaction area underneath are thoroughly wet by urine, you should replace it with a new test. If the area is still dry, clean away the urine with tissue paper. Follow the instructions to run the test. If a purple band does not appear in the C area in 5 minutes, the test should be voided. However, if you do get a band in the control window, you can rely on your results.

* 99% accurate at detecting typical pregnancy hormone levels. However, some pregnant women may not have detectable amounts of pregnancy hormone in their urine on the first day of the missed period or may have miscalculated the first day of their period.

Do not use a test that is past its expiration day. The test is not reusable. Do not use if the foil pouch is damaged. Store in a dry place 36-86 °F (2-30°C). Do not freeze. For In-Vitro diagnostic use only. Not for internal use. Keep out of the reach of children. Dispose of properly after use.

